





1. Dashboard

- Overview of today's cases (assigned, pending, reported)
- Notifications/alerts (urgent studies, revisions requested)
- Quick access to unread studies
- Filters by name, modality, date, referring doctor, etc.

2. My Cases

- Patient demographics (ID, name, age, sex)
- Study details (modality, date/time, body part, reason for exam)
- Status (pending, in review, reported)
- Priority tagging
- Advanced filtering and search
- Sort by: Name, modality, referring physician, urgency, etc.

The screenshot displays the 'My Cases' interface in the MEDICALAI application. The top navigation bar includes the MEDICALAI logo, a dropdown menu for 'Medi Care - De...', and an 'Invite to Medical' button. A notification bell icon with a red '1' and a user profile for 'Shan Marken' are also present. The main content area is titled 'My Cases' and shows a search bar with 'smith' entered. Below the search bar are filters for 'Labels', 'Case Manager', 'Members', and 'Institution name'. A 'Last modified' date filter is also visible. A table lists two cases, both created on 17/06/2025. The first case is for 'smith' and has a status of '+1'. The second case is also for 'smith' and has a status of '+1'. A notification at the bottom of the table states: 'Last update: Dr. Michael Jude invited you to participate in the case of smith.' The interface includes a sidebar with various icons for navigation and a 'New case +' button in the top right corner.

Patient details	Patient ID	Members	Created at	Actions
<input type="checkbox"/>  smith	-	 +1	17/06/2025	Share More
Last update: Dr. Michael Jude invited you to participate in the case of smith.				
<input type="checkbox"/>  smith	-	 +1	17/06/2025	Share More

Patient Studies

+ New Study

ID	Name	Report	Accession	Modality	Description	Date	Time	Source AE	Actions
PID001	John Doe	Report_001.pdf	ACC-203948	CT	Head CT Scan	2024-12-01	10:15AM	CT-SOURCE-AE	
PID001	John Doe	Report_001.pdf	ACC-203948	CT	Head CT Scan	2024-12-01	10:15AM	CT-SOURCE-AE	
PID001	John Doe	Report_001.pdf	ACC-203948	CT	Head CT Scan	2024-12-01	10:15AM	CT-SOURCE-AE	
PID001	John Doe	Report_001.pdf	ACC-203948	CT	Head CT Scan	2024-12-01	10:15AM	CT-SOURCE-AE	
PID001	John Doe	Report_001.pdf	ACC-203948	CT	Head CT Scan	2024-12-01	10:15AM	CT-SOURCE-AE	
PID001	John Doe	Report_001.pdf	ACC-203948	CT	Head CT Scan	2024-12-01	10:15AM	CT-SOURCE-AE	
PID001	John Doe	Report_001.pdf	ACC-203948	CT	Head CT Scan	2024-12-01	10:15AM	CT-SOURCE-AE	

10 per page

Showing 1-10 of 77 studies

< Previous

1

2

3


...

8

Next >

3. DICOM Viewer Integration

- Seamless link to the DICOM viewer




4. My workspace

All the pending, in review cases will be appeared here


Patient demographics (ID, name, age, sex)

Study details (modality, date/time, body part, reason for exam)




Reporting Module


- Structured and free-text reporting templates



- Template management (create/save)






- **Insert key images into reports**
- Report auto-save and versioning
- Auto population of patient & study metadata



- Reporting Workflow

- Draft / 2nd opinion / Finalize / Sign-off options
- Second opinion request (Image & report sharing via secure link)
- Internal notes/comments per case
- Add signature





2. Admin (User Management & Permissions)

- Roles: Radiologist
- Create new Cases
- Give Access control per study
- Ability to delegate or reassign cases



3. Dentist

- Patient demographics (ID, name, age, sex)
- Study details (modality, date/time, body part, reason for exam)
- View Image and Report

1D

3D

1W

1M

1Y

ALL

Search...

Reset

Search

Modality :

AI

CT

MRI

DX

IO

CR

Source : All Sources

Patient Studies

10 per page

Showing 1-10 of 77 studies

< Previous

1

2

3

...

8

Next >

4. Patient

- Patient demographics (ID, Name)
- Study details (Description, Modality, Date, Time, Action)
- View Image and Report

Pages / Home

CARE FOR YOUR SMILE

MAKE AN APPOINTMENT NOW

Scans and X rays

Date : 2023-10-23
Time: 10.23 A.M

Date : 2023-10-23
Time: 10.23 A.M

Lab Reports

Date : 2023-10-23
Time: 10.23 A.M

Date : 2023-10-23
Time: 10.23 A.M

Date : 2023-10-23
Time: 10.23 A.M

Stories